



PRECISION[™]
SURGICAL WORKS

DENTAL ORTHODONTIC INSTRUMENTS


Zahnärztliche Diagnostikinstrumente
Instrumentos de Dentales Diagnóstico
Instruments de Dentaire Diagnostic
Dentale Diagnostica Strumenti



Multi Use Cutter

For Soft Wire Pins and Ligature Wires upto 0.38 mm

Catalog #	Size
PRO-32-989-01	Standard



Mini Pin and Ligature Cutter

For Soft Wire Pins and Ligature Wires upto 0.30 mm

Catalog #	Size
PRO-32-989-02	Standard



Micro Mini Pin and Ligature Cutter

For Soft Wire Pins and Ligature Wires upto 0.30 mm

Catalog #	Size
PRO-32-989-03	Standard



Micro Mini Pin and Ligature Cutter

With Long Handle - For Soft Wire Pins and Ligature Wires upto 0.30 mm


Catalog #	Size
PRO-32-989-04	Standard



Micro Mini Pin and Ligature Cutter

Angled 15° - For Soft Wire Pins and Ligature Wires upto 0.30 mm

Catalog #	Size
PRO-32-990-05	Standard



Universal Distal End Cutter

With Long Handle - For Cutting All Types of Wires from 0.30 mm to 0.53 mm x 0.64 mm

Catalog #	Size
PRO-32-990-01	Standard



Universal Distal End Cutter

With Long Handle - For Cutting All Types of Wires from 0.30 mm to 0.53 mm x 0.64 mm

Catalog #	Size
PRO-32-990-02	Standard





Flush Cut Distal End Cutter No Hold
For All Types of Wires up to a maximum of 0.53 mm x 0.64 mm

Catalog #	Size
PRO-32-990-03	Standard




Wire Cutter
Angled 15° - For Hard Wires up to 0.56 mm x 0.71 mm

Catalog #	Size
PRO-32-990-04	Standard



Wire Cutter
Straight - For hard Wires up to 0.56 mm x 0.71 mm

Catalog #	Size
PRO-32-990-05	Standard



Wire Cutter

For Cutting Soft Wire Pins and Ligature Wires up to 0.30 mm

Catalog #	Size
PRO-32-990-06	Standard



Flush Cut and Hold Distal End Cutter

For Cutting All Types of Wires up to 0.51 mm

Catalog #	Size
PRO-32-990-07	Standard



Weingart Utility Plier

Catalog #	Size
PRO-32-990-08	Standard



**Weingart Utility Plier
Slim**

Catalog #	Size
PRO-32-990-09	Standard

**How Utility Plier
Straight**

Catalog #	Size
PRO-32-990-10	Standard



How Utility Plier
Curved

Catalog #	Size
PRO-32-990-11	Standard



Band Seating Plier

Catalog #	Size
PRO-32-990-12	Standard



Adhesive Removing Plier

Catalog #	Size
PRO-32-990-13	Standard



Posterior Band Removing Plier

Catalog #	Size
PRO-32-990-14	Standard



Posterior Band Removing Plier

Catalog #	Size
PRO-32-990-15	Standard



Tip Back Plier

Catalog #	Size
PRO-32-990-16	Standard



Force Module Separating Plier

Catalog #	Size
PRO-32-990-17	Standard



Coon Suture Tying Plier

Catalog #	Size
PRO-32-990-18	Standard



Crown & Band Contouring Plier

Catalog #	Size
PRO-32-990-19	Standard



Bracket Removing Plier

Catalog #	Size
PRO-32-990-20	Standard





Bracket Removing Plier
Angulated

Catalog #	Size
PRO-32-990-21	Standard



Band Removing Plier

Catalog #	Size
PRO-32-990-22	Standard



Jarabak Wire Forming Plier

Catalog #	Size
PRO-32-990-23	Standard



Three Jaw Plier

Catalog #	Size
PRO-32-990-24	Standard



Ortho Finishing Step Plier

Catalog #	Size
PRO-32-990-25	Standard



Omega Loop Forming Plier

Catalog #	Size
PRO-32-990-26	Standard

Bird Beak Plier

Catalog #	Size
PRO-32-990-27	Standard



Light Wire Bird Beak Plier

Catalog #	Size
PRO-32-990-28	Standard



Tweed Loop Forming Plier

Catalog #	Size
PRO-32-990-29	Standard



Arch Bending Plier

Catalog #	Size
PRO-32-990-30	Standard



**Tweed Plier
Slim**

Catalog #	Size
PRO-32-990-31	Standard



Rectangular Arch Bending Plier

Catalog #	Size
PRO-32-990-32	Standard



Lingual Arch Forming Plier

Catalog #	Size
PRO-32-990-33	Standard



Hollow Chop Contouring Plier

Catalog #	Size
PRO-32-990-34	Standard



O'Brien Plier

Catalog #	Size
PRO-32-990-35	Standard



Mini Three Jaw Plier

Catalog #	Size
PRO-32-990-36	Standard



Marcotte Looping Plier

Catalog #	Size
PRO-32-990-37	Standard



Hollow Chop Plier

Catalog #	Size
PRO-32-990-38	Standard



Light Wire Bird Beak Plier

Catalog #	Size
PRO-32-990-39	Standard



Nance Loop Forming Plier

Catalog #	Size
PRO-32-990-40	Standard



Adams Plier

Catalog #	Size
PRO-32-990-41	Standard



V-Stop Plier

Catalog #	Size
PRO-32-990-42	Standard



Hook Crimping Plier

Catalog #	Size
PRO-32-990-43	Standard



Optical Plier

Catalog #	Size
PRO-32-990-44	Standard



Loop Forming Plier

Catalog #	Size
PRO-32-990-45	Standard



Bird Beak Plier With Cutter

Catalog #	Size
PRO-32-990-46	Standard



Long Tapered Bird Beak Plier

Catalog #	Size
PRO-32-990-47	Standard



Hammerhead Tie Back Plier

Catalog #	Size
PRO-32-990-48	Standard





Ligature Director

Catalog #	Size
PRO-32-993-01	Standard



Catalog #	Size
PRO-32-993-02	Standard



Ligature Director/Pluggger

Catalog #	Size
PRO-32-993-03	Standard



Tucker Ligature Director

Catalog #	Size
PRO-32-993-04	Standard



Plugger/Ligature Picker

Catalog #	Size
PRO-32-993-05	Standard



Fillion Explorer/Ligature Director

Catalog #	Size
PRO-32-993-06	Standard



Mershon Band Pusher

Catalog #	Size
PRO-32-995-01	Standard



Band Seater

Catalog #	Size
PRO-32-995-02	Standard



Band Pusher / Scalar
Fig. 1

Catalog #	Size
PRO-32-997-01	Standard



Band Pusher / Scalar
Fig. 2

Catalog #	Size
PRO-32-997-02	Standard



Direct bond Bracket Holder

Catalog #	Size
PRO-32-998-00	Standard